

## Prosthetic Materials Needed for Day of Surgery

In preparation for the delivery of an nSequence® Immediate Load Long-Term Provisional prosthesis, it is strongly recommended that you have the following restorative materials available on the day of the surgical procedure. These materials are not provided by nSequence:

- **Burs** – A small to medium size “football shaped” cross-cut bur for smoothing the provisional as needed and adjusting the receptor sites on the long term provisional). A #4 Round Bur, #1557 Carbide Bur (or similar) and a #1 Round Bur for a slow speed handpiece, are also helpful for refining small areas and opening any air bubbles.
- **Separating Disc** – For quick and easy removal of the arms on the Long-Term Provisional. Suggest; Komet®, Honeycomb Design, 6934B.
- **Bone Reduction Bur** – Suggested size/shape is a; Komet, S-Cutters H72E HP, H72E.104.060, or similar.
- **Rubber Points (2-3 each)** – Medium coarse and fine coarse rubber points (for smoothing acrylic prior to polishing). If you would like to purchase a polishing kit that does not require pumice, we would recommend the Kerr® Acrylic Adjustment Polishing Set part number LS-901. It includes “football shaped” acrylic burs, as well.
- **Luting/Bonding Material** – We would recommend Quick Up™ by Voco® or Chairside® from Zest® Dental Solutions.
- **Teflon® Tape and Composite Material** – Teflon Tape to plug the screw access holes prior to sealing these access holes with Fermit® by Ivoclar Vivadent® composite.
- **Handpiece/Motor** – For use with the Laboratory Burs (Straight Burs not Latched Burs).
- **Can of Compressed Air** – Same as used for cleaning off computer keypads.
- **Soft Reline Material** – Used to help ensure a secure fit of the backup denture, in situations when the long-term provisional cannot be delivered on the day of surgery.
- **Horseshoe Articulating Paper** – For initial occlusal adjustments.
- **Light Cure Unit** – Hand-held for dual cure resins.
- **Petroleum Jelly** – To act as a release agent for self-curing resin.

**For more information, please contact nSequence at 888.809.2777  
or visit us at [www.nsequence.com](http://www.nsequence.com).**

\*These are recommended materials. Ultimately, each implant team should use their own discretion regarding what materials they choose to have available on the day of surgery.