

NEW INNOVATIONS = NEW **TIME-SAVING** PROTOCOLS

Say goodbye to Prosthetic Bite Guides (PBG) and get faster results with our newly-designed Long Term Provisionals.



IMPROVED TECHNIQUE

Open Mouth Luting without a PBG

- ✓ Easier to lute sedated patients
- ✓ Easier and faster to light cure material
- ✓ Precisely shape/contour occlusal with a perio probe
- ✓ No need to hold your patient in occlusion

Please contact us for more information at [nSequence.com](https://www.nSequence.com) or **888.809.2777**

Open Mouth Luting without a Prosthetic Bite Guide



One injection hole per site



Arms secure Provisional to
Bone Foundation Guide



Full occlusal view for luting



Fill from apical up to occlusal



Perio probe used to contour occlusal



Results in less finishing time

Please contact us for more information at [nSequence.com](https://www.nSequence.com) or 888.809.2777