

TISSUE SUPPORTED		BONE SUPPORTED	
WITH PROVISIONAL	WITHOUT PROVISIONAL	WITH PROVISIONAL	WITHOUT PROVISIONAL
FULLY DENTATE		FULLY DENTATE	
<ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Upper and Lower Impressions ✓ Bite Registration 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Working Arch Impression
PARTIALLY DENTATE		PARTIALLY DENTATE	
<ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Physical Wax Rims ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Physical Wax Rims ✓ Upper and Lower Impressions ✓ Bite Registration 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Physical Wax Rims ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Physical Wax Rims ✓ Working Arch Impression
FULLY EDENTULOUS		FULLY EDENTULOUS	
<ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Denture or Wax Rim ✓ Send Physical Wax Rims (if applicable) ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Denture or Wax Rim ✓ Send Physical Wax Rims (if applicable) ✓ Upper and Lower Impressions ✓ Bite Registration 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Denture or Wax Rim ✓ Send Physical Wax Rims (if applicable) ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Denture or Wax Rim ✓ Send Physical Wax Rims (if applicable) ✓ Working Arch Impression

TEETH AND TISSUE		TEETH AND BONE	
WITH PROVISIONAL	WITHOUT PROVISIONAL	WITH PROVISIONAL	WITHOUT PROVISIONAL
FULLY DENTATE		FULLY DENTATE	
<ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Working Arch Impression 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Working Arch Impression
PARTIALLY DENTATE		PARTIALLY DENTATE	
<ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Wax Rim with Markers ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Wax Rim with Markers ✓ Working Arch Impressions 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Wax Rim with Markers ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos 	<ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Wax Rim with Markers ✓ Working Arch Impressions