

| TISSUE SUPPORTED | | BONE SUPPORTED | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| WITH PROVISIONAL | WITHOUT PROVISIONAL | WITH PROVISIONAL | WITHOUT PROVISIONAL |
| FULLY DENTATE | | FULLY DENTATE | |
| <ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Upper and Lower Impressions ✓ Bite Registration | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Working Arch Impression |
| PARTIALLY DENTATE | | PARTIALLY DENTATE | |
| <ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Physical Wax Rims ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Physical Wax Rims ✓ Upper and Lower Impressions ✓ Bite Registration | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Physical Wax Rims ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Physical Wax Rims ✓ Working Arch Impression |
| FULLY EDENTULOUS | | FULLY EDENTULOUS | |
| <ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Denture or Wax Rim ✓ Send Physical Wax Rims (if applicable) ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Denture or Wax Rim ✓ Send Physical Wax Rims (if applicable) ✓ Upper and Lower Impressions ✓ Bite Registration | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Denture or Wax Rim ✓ Send Physical Wax Rims (if applicable) ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Denture or Wax Rim ✓ Send Physical Wax Rims (if applicable) ✓ Working Arch Impression |

| TEETH AND TISSUE | | TEETH AND BONE | |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| WITH PROVISIONAL | WITHOUT PROVISIONAL | WITH PROVISIONAL | WITHOUT PROVISIONAL |
| FULLY DENTATE | | FULLY DENTATE | |
| <ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Upper and Lower Impressions ✓ Bite Registration | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient ✓ Working Arch Impression |
| PARTIALLY DENTATE | | PARTIALLY DENTATE | |
| <ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Wax Rim with Markers ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Wax Rim with Markers ✓ Upper and Lower Impressions ✓ Bite Registration | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Wax Rim with Markers ✓ Upper and Lower Impressions ✓ Bite Registration ✓ Photos | <ul style="list-style-type: none"> ✓ CBCT Scan of Patient Wearing Wax Rim with Markers ✓ Send Wax Rim with Markers ✓ Working Arch Impressions |